

User Manual for

Automated Medicine Dispenser

“Your Daily Health Partner”

MedSmart Innovations

31st October 2023

Welcome to the user manual for MedSmart Innovations' Medicine Dispenser.

We are delighted to introduce you to our innovative solution designed to make managing your medication a breeze. This manual is your guide to using and maintaining your medicine dispenser effectively. It contains step-by-step instructions, helpful tips, and all the information you need to ensure a seamless and safe experience. We understand the importance of your health, and our product is here to simplify the process of medication management for you.

Let's get started!



Contents

Safety Precautions:.....4

Getting Started and Using the Medicine Dispenser:.....5

Setting the Time with Push Buttons:.....6

Technical Specifications:.....7

Contact Information:.....7

Safety Precautions:

Keep Out of Reach: Ensure that the medicine dispenser is kept out of reach of children and pets. It contains medication that must not be consumed inadvertently.

Use as Directed: Only use the medicine dispenser as instructed in this manual. Follow your healthcare provider's medication regimen precisely.

Storage: Store the medicine dispenser in a clean, dry place away from extreme temperatures, moisture, and direct sunlight.

Repairs: Do not attempt to repair the dispenser yourself. Contact our customer support or a qualified technician for any necessary repairs.

Emergency Situations: In case of emergencies or if you suspect a malfunction, discontinue use and seek immediate medical advice.


Medication Labeling: Ensure that medications are correctly labeled, and the dispenser is loaded with the right medicines. Double-check before using.

Children and Seniors: If the dispenser is intended for use by children or seniors, consider their specific needs and abilities. Simplify usage where necessary.

Medication Allergies: Be aware of any allergies you or the user may have to medications. Check each medication carefully.

Getting Started and Using the Medicine Dispenser:

Open the Top Cover: Lift the top cover of the medicine dispenser.



Load Medication: Inside, you'll find eight compartments, each corresponding to the scheduled doses. Place the required medication for each dose into these compartments.

Close the Top Cover: Carefully close the top cover as you did before, ensuring it latches securely.

Power Up: Connect the medicine dispenser to a power source using the provided power cord.

Set the Real-Time Clock: Use the push buttons located below the LCD screen to set the current time on the device.

Choose Time or Alarm Mode: The Time/Alarm switch allows you to select between setting the real-time clock or configuring alarms.

Set Alarms: If you've selected the alarm mode, use the buttons to set the specific times for medication reminders.

When the Alarm Sounds: When an alarm triggers, the medicine dispenser will release the scheduled medication into the designated drawer. The patient can then easily access and take the medicine from the drawer.

Setting the Time with Push Buttons:

- To accurately set the time, the medicine dispenser is equipped with four push buttons. These buttons allow you to increase or decrease the hours and minutes on the device's real-time clock. The buttons are typically labeled as "+ Hr" and "- Hr" for hours, and "+ Min" and "- Min" for minutes.
- When you press the "+ Hr" button, it will increase the hour value. Conversely, the "- Hr" button will decrease it. Likewise, the "+ Min" and "- Min" buttons serve to adjust the minutes.
- Make sure to use these buttons with precision when setting the real-time clock and configuring alarms. They provide an intuitive way to ensure that the medicine dispenser accurately reflects your schedule.

Technical Specifications:

Dimensions:

Height: 26 cm

Width: 12 cm

Weight: 318g

Power Input: 5V

Capacity: 7 doses

Display: LCD 16x2 display

Contact Information:

If you have a question about the device or to report any errors:

Contact us:

077-8797936

2020s18114@stu.cmb.ac.lk